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
Figure 3: Tubestand (Gantry with C-arm) Measurements

A.	Height	223 cm (87.8 inches)
B.	Width	66 cm (26 inches)
C.	Width	173 cm (68.0 inches)
D.	Depth	138 cm (54.3 inches)
	Weight	Maximum of 400 kg (882 pounds)

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
Figure 5: Premium Acquisition Workstation Measurements

A.	Height	202 cm (79.8 inches)
B.	Width	92.7 cm (36.5 inches)
C.	Depth	58.5 cm (23.0 inches)
	Weight	154 kg (340 pounds)

Room Size

Length	365 cm (12 feet)
Width	275 cm (9 feet)
Ceiling	244 cm (8 feet)
Door	214 cm x 92 cm (84 in. x 36 in.)

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Notes:

The exam room layout should be pre-planned before the arrival of the system. Refer to [Room Layout Worksheet](#) on page 44. A typical room layout is shown here. Check specific local or hospital requirements for additional data. Check the route from the loading dock to the room.

Maintain a clearance of 50 cm between the tube head at any position of rotation and any stationary object (excluding the floor). If this clearance cannot be maintained due to space constraints, provide as much clearance as possible. The standard operation of the controls for C-arm rotation should not be altered at installations that do not meet the 50 cm clearance.


Figure Legend

1. Universal Acquisition Workstation
2. Tubestand/Gantry
3. Minimum clearance of 50 cm (20 in.)
4. Consult local regulations for minimum clearance
5. Minimum doorway opening of 92 cm (36 in.)
6. Patient Area

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Improving the early detection and diagnosis of **breast cancer**.

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
Image is for illustrative purposes. Actual software solution view will be different.

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Integrated Artificial Intelligence Solutions

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Genius AI® Detection Solution

What if you could detect almost one additional cancer?^{1**}

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“

The 3DQuorum algorithm lessens what would otherwise be a voluminous and untenable dataset into one more palatable for our PACS and network systems and **more easily interpreted by radiologists with reduced turnaround time, fatigue, and digital eye strain.**”

Jason McKellop, MD

Director of Women's Imaging,
Beverly Radiology Medical Group and
Medical Director, Breastlink Tarzana

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Contrast Enhanced Mammography (CEM)
powered by I-View 2.0TM Software

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Contrast Enhanced Mammography (CEM)

Contrast-enhanced mammography (CEM) uses IV injection of a dye containing iodine (as is used in CT) in combination with a standard digital mammogram.


CEM is FDA-approved for diagnostic breast imaging, and is proven effective in looking at the extent of cancer in women with newly diagnosed breast cancer, to check a cancer's response to chemotherapy before surgery, and to assess breast symptoms, or abnormalities seen on standard mammograms.

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Contrast Enhanced Mammography (CEM)

CEM vs. MRI

	Faster Screening Time	Higher Patient Compliance	Better Patient Experience
MRI	30-60 min. imaging 60-90 min. room time ⁵	Only 58% compliance ⁶ , time-consuming for pre-authorization insurance process	Anxiety-provoking, long-time laying prone, discomfort ⁷
CEM	8-20 min. imaging time 30 min. room blocking time	25% of the cost of MRI ⁵ , radiologists can order immediately without pre-authorization	79% of patients prefer CEM over MRI ⁷ , same positioning as standard mammogram

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Contrast Enhanced Mammography (CEM)

Standard

- I-View™ 2.0 CEM Imaging Software License (single gantry)
- Copper filter

Requirements

- Selenia® Dimensions® or 3Dimensions™ with 1.10/2.1 software or higher
- Diagnostic license

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Contrast Enhanced Mammography (CEM)

Site preparation


Your medical physicist must perform the necessary testing for CEM

Obtain a portable power injector with contrast material

A chair for the patient – comfort for administration of contrast

- For the imaging portion, the patient may stand or sit

Protocol for the procedure is determined by the Radiologist

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Contrast Enhanced Mammography (CEM)

Quality Control Tests for I-View™ procedure

One additional Gain Calibration
(Detector Flat Field Calibration)

- Required weekly
- CEM Gain Calibration

One additional view during the Artifact
Evaluation

- Required weekly
- For the copper filter

Thank you for your attention

If you need any further
information, Please contact us

Tel: +982188811106